**DESIONIZADOR DE AGUA LAB-ION L2
DEMINERALIZER OF WATER LAB-ION L2
DEMINERALISATEUR D'EAU LAB-ION L2**

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ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish.....	2-3
English.....	4-5
French	6-7

1. Remove cartridge and conductivity meter from the brackets, then mount brackets onto the wall near your water supply, using the screws and expanding plugs provided.

2. Place cartridge and conductivity meter back into the brackets. Connect water inlet hose by applying slight pressure against the cartridge, followed by a quarter turn of the quick coupling until it rests firm place.

Please, take care that the inlet hose has no twisted and runs straight from the tap to the inlet connection.

3. Now push the water outlet hose over the nozzle on the base of the conductivity meter.

4. Plug into your power supply with the plug provided (if a different type of power plug is used in your country, please modify). The unit is now ready to produce de-ionised water whenever you turn your tap on.

NOTE: Do not alter the bore of the reduction adapter which fits onto the water tap.

INSTRUCTIONS FOR EXCHANGE OF CARTRIDGES

By watching the readings on your conductivity meter you can determine when the resin is exhausted and the cartridge needs to be exchanged.

Please, follow these simple instructions:

1. Turn off water tap, then disconnect from power supply. Disconnect water inlet hose from the water tap and empty the unit of water by placing the inlet hose into a lower situated container of basin.

2. Remove the unit from the brackets, disconnect the exhausted cartridge from the conductivity meter and dispose. The contents of the exhausted cartridge may be used as a fertilizer or aerat for soil.

3. Reconnect the conductivity meter onto the new cartridge and place back into the brackets, then re-connect the inlet and outlet hoses as described above in the installation guide, paragraph 2 and 3.

MAINTENANCE

1. After a certain time of usage of the LAB-ION it is possible that the sieve (which is located under the conductivity check) is plugged from residues of tap water or little parts of resin. We recommend in turning out and cleaning the sieve under running water, when you will change the exhausted cartridge.

2. If you do not check the sieve at times it is possible that the sieve will be totally plugged and the demineralized water has no free run. In this case high pressure in the cartridge can damage the unit (hoses can get loose...).

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.